

Dr. Brown's Natural Flow® Baby Bottle

CLEANING


The cleaning of the two (2) additional vent assembly parts will require an extra moment. This is time well spent to help avoid spit-ups and other discomforts to your baby.

1. Always use clean hands and surfaces.
2. Fully disassemble all parts.
3. **Silicone Nipple - first use only:** boil the nipple in water for three (3) minutes.
4. Wash all parts thoroughly in hot water and detergent. Use enclosed brush on all holes in the insert and reservoir.
5. Rinse parts thoroughly, making sure that the nipple hole is clean and clear.
6. If you use a dishwasher, keep all parts in the top rack only. Do not exceed 140°F (60°C) water temperature.

ASSEMBLY and FEEDING INSTRUCTIONS (Please refer to diagram)

NOTE: The bottle WILL NOT FUNCTION without the INSERT and RESERVOIR.

1. **Avoid formula bubbles:** Do not shake the bottle. Formula should be stirred gently to fully dissolve lumps - they can clog nipple holes.
2. Fill to the desired amount. The vent will not work properly if you fill the bottle above the FILL LINE WARNING.
3. Snap the reservoir fully onto the insert.
4. Place the reservoir into the bottle, with the insert making full contact with the top of the bottle.
5. Snap the nipple into the nipple collar and place loosely on the bottle.
6. Warm the bottle - Always remove the nipple and collar before warming to prevent the vent from leaking fluid if gas forms during warming.
7. Do not warm the liquid beyond body temperature of 98.6°F (37°C).
8. **Feed - Be sure to ALWAYS tighten the collar firmly before feeding.**



CAUTIONS

1. Due to the risk of creating formula hot spots and overheating, the use of a microwave oven is not recommended. If you use a microwave oven, remember that the nipple and collar must always be removed during heating. Test the temperature frequently by placing a few drops on your tongue or back of your hand after stirring thoroughly to eliminate hot and cold spots.
2. Pull on the nipple before each use to make sure it does not tear. Replace the nipple if it shows any signs of wear, such as cracking or change of shape.
3. Stop feeding before your baby has completely emptied the nipple to keep your baby from swallowing air that may cause colic or other digestive problems.
4. Remove the bottle from your baby's mouth when the baby is not feeding. Prolonged contact with milk, formula, and other liquids can cause tooth decay.
5. The components of this bottle are not toys and should not be given to children to play with or suck on.
6. Be sure to contact your physician should you have any feeding problems or questions.

HANDI-CRAFT COMPANY
4433 Fyler, St. Louis, MO 63116
U.S.A.

To order bottles and nipples call:
1-800-778-9001

For more information visit our website at:
www.handi-craft.com

The product is protected by two or more U.S. patents (5,570,796 5,779,071 D420,448 D423,107) and/or national patents, with other U.S. and national patents pending.

REMINDERS

A. The nipple included has been specifically selected for an infant's first use of the free-flowing characteristics of the Dr. Brown's Natural Flow® venting system. Faster flowing Stage 2 and Y-Cut/Cereal versions of the nipple are available. If Dr. Brown's™ nipples are not available from your retailer, they can be ordered directly from Handi-Craft Company. Nipples of other brands may leak or feed too fast or too slow.

B. When feeding your baby, hold the bottle so that the reservoir tube hole in the bottom of the bottle is not covered by liquid. (Approximately at a 45° angle)

C. Keep your baby propped up rather than lying down during feeding.

D. While the vent continuously equalizes pressure, movement of any liquid remaining in the reservoir tube is normal.

E. If the reservoir bulb remains over 1/3 full when the bottle is placed upright, please check that the nipple hole is clear and that the insert and reservoir are firmly pressed together. The reservoir bulb can be emptied by momentarily loosening the nipple collar. (Be sure to retighten before feeding)

F. When transporting bottles, always keep them in the upright position.

G. For proper bottle and venting functions use only Dr. Brown's™ nipples.

02-04-30
©Handi-Craft Company 2002